


HAYDEN MEDICAL

User Manual

FULL HD Endoscope System Manual



1. Functions of Keys and Interfaces



1. 19" HD Medical Monitor

2. Screen Operation Panel

3. Display ON/OFF

4. Settings Panel

5. Light Source Interface

6. Camera Head Interface

7. USB Storage

8. Power ON/OFF

9. COM 2

10. JMICK1

11. COM1

12. SIM-1

13. LED POWER

14. USB

15. CVBS Video Output

16. DVI Video Output

17. HDMI Video Output

18. Grounding Rod

19. Power Switch

20. Fuse

21. Input, (AC 110V-220V)


22. Image Capture

23. Freeze

24. Video Recording

25. White Balance

2. Functions of Operating Panel




1. Up

2. Down

3. Left (useless)

4. Right (Adjust monitor color pattern)

5. Menu/Confirmation



1. Lock Menu

2. Left for Controlling Menu

3. Image Capture

4. Up for Controlling Menu

5. Down for Controlling Menu

6. Camera Menu

7. White Balance

8. Right for Controlling Menu

9. Sleep of Light Source

10. Video Recording

11. Reduce LED Light

12. Increase LED light

13. Switch Over to Software

14. LED Light Source Indicator Light

15. Machine usage time

3. Installation for Reference





4. Video Recording



(1) Image Capture

1. Insert the USB Flash Drive into the USB Port on the front of the machine
2. Press **[F] 22** button (Image Capture Button) on camera head to capture on-screen image
3. The **[Capture Indicator Light]** will illuminate on the operation panel
4. The image is saved onto the USB flash drive and **[Capture Indicator Light]** turns off

(2) Video Recording

1. Insert the USB Flash Drive into the USB Port on the front of the machine
2. Press **[Recording]** button on operation panel to start video capture
3. The **[Capture Indicator Light]** will illuminate next to the USB Port
4. The **[Capture Indicator Light]** will remain lit when video is being captured
5. Press the **[Recording]** button again to end video capture
6. The **[Capture Indicator Light]** flickers 3 times while the video is being saved
7. The video is saved onto the USB Flash Drive, and the **[Capture Indicator Light]** turns off