

MOBILE SURGI-STRETCHER



HYDRAULIC UNITS FOR OPHTHALMOLOGY

Distributed By:

HAYDEN MEDICAL

661-296-7200

info@haydenmedical.com

www.haydenmedical.com

www.haydenophthalmic.com



578 SURGI-STRETCHER SERIES

The Surgi-Stretcher Series is designed as a mobile surgical platform for EYE/ENT and plastic surgery procedures.



VERSATILITY:

- Articulating head piece with dual operating knobs allow infinite adjustability and precise microscopic movement of the head section
- Designed with a low-profile head piece allowing site access for both superior (over brow) and lateral (side approach) procedures



STANDARD FEATURES:

- Four wheel brake and steer caster system operable from all four corners
- Easy to clean thermoformed base cover with O2 tank holder and ample storage space
- 3-Stage Telescoping I.V. pole [permanently mounted at foot]
- Foot-end hand crank controlled backrest and leg flex
- Stainless steel Retracto[®] Side Rails with rail caps
- Dual side hydraulic foot pedal controls
 - Height Adjustment
 - Trendelenburg Positioning
- Comfort mattress and head pad
- Dual operating control knobs
- Easy rolling casters
- Low profile articulating head piece



Microscopic articulation of the head piece is achieved by rotating the control knobs to raise/lower the patient's forehead or chin.



Specifications:

	578EYEST 578EYXST	
	inch	cm
Overall Width:	29.75	75.6
Overall Width: <i>EYX models only</i>	34.5	87.6
Overall Length:	88.5	224.8
Overall Height: <i>range - minimum</i>	24	61
Overall Height: <i>range - maximum</i>	32	81.3
<i>patient surface</i> Width:	25	63.5
<i>patient surface</i> Width: <i>EYX models only</i>	30	76.2
<i>patient surface</i> Length:	81.5	207
<i>retracto</i> [®] rail Length:	49	124.5
<i>retracto</i> [®] rail Height:	14	35.6
<i>comfort</i> Mattress Pad:	3	7.6
<i>easy rolling</i> Casters:	8	20.3
	degree of movement: ²	
Backrest: <i>movement</i>	0 - 90	
Knee Flex: <i>movement</i>	0 - 75	
Trendelenburg: <i>standard</i>	16	
Trendelenburg: <i>reverse</i>	16	
Weight Capacity: <i>EVENLY DISTRIBUTED</i>	lb	kg [†]
	500	227
Shipping Weight: <i>EYEST</i>	410	186
Shipping Weight: <i>EYXST</i>	450	204

Features:

Operation:	Hydraulic
<i>dual side hydraulic</i> Foot Pedal Controls: ²	■
<i>foot-end hand crank controlled</i> Backrest:	■
<i>foot-end hand crank controlled</i> Leg Flex:	■

■ Standard — Not Applicable/Available

² foot pedal controls adjust height & trendelenburg positioning

[†] standard rounding method used to show metric values as whole numbers

578 Product Info



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

Distributed By:

HAYDEN MEDICAL

661-296-7200

info@haydenmedical.com

www.haydenmedical.com

www.haydenophthalmic.com

Product Hyperlinks Enabled

graham field

GF1500008RevE22