

MOBILE SURGI-STRETCHER



ELECTRIC UNITS FOR OPHTHALMOLOGY

Distributed By:
HAYDEN MEDICAL
661-296-7200
info@haydenmedical.com
www.haydenmedical.com
www.haydenophthalmic.com

5E82 ELECTRIC SURGI-STRETCHER SERIES

The Surgi-Stretcher Series is designed as a mobile surgical platform for EYE/ENT and plastic surgery procedures.



VERSATILITY:

- Articulating head piece with dual operating knobs allow infinite adjustability and precise microscopic movement of the head section
- Designed with a low-profile head piece allowing site access for both superior (over brow) and lateral (side approach) procedures
- Battery powered stretcher with easy to use interactive hand pendant.
- Leg flex operated controls are available for height adjustment, backrest, leg section, and Trendelenburg/Reverse Trendelenburg



Specifications:

	5E82EYEST <small>standard width</small>		5E82EYXST <small>wide width</small>		
	inch	cm†	inch	cm†	
Overall Width:	32.00	81.28	37.00	93.98	
Overall Length:	92.34	234.54	92.34	234.54	
Overall Height: <small>minimum</small>	22.00	55.88	22.00	55.88	
Overall Height: <small>maximum</small>	40.00	101.60	40.00	101.60	
Patient Surface: <small>width</small>	25.00	63.50	30.00	76.20	
Patient Surface: <small>length</small>	85.38	216.87	85.38	216.87	
Side Rail Length:	48.00	121.92	48.00	121.92	
Side Rail Height:	14.00	35.56	14.00	35.56	
Mattress Pad: <small>thickness</small>	3.00	7.62	3.00	8.762	
Casters:	5.90	14.99	5.90	14.99	
		degrees °		degrees °	
Backrest: <small>movement</small>		0 - 85		0 - 85	
Knee Flex: <small>movement</small>		0 - 60		0 - 60	
Trendelenburg: <small>standard</small>		15		15	
Trendelenburg: <small>reverse</small>		15		15	
		power		power	
Volts (V)		100 - 240		100 - 240	
Hertz (Hz)		~ 50/60		~ 50/60	
		lb	kg†	lb	kg†
Weight Capacity: ¹ <small>EVENLY DISTRIBUTED</small>		800	362.87	800	362.87
Shipping Weight:		420	190.51	460	208.65

STANDARD FEATURES:

- Four wheel brake and steer with electric casters system
- Base cover with built in O2 tank holder and ample storage space
- 3-Stage Telescoping IV pole [permanently mounted at foot]
- Stainless steel Retracto® Side Rails with rail caps
- Low-profile articulating head piece with dual operating control knobs
- Comfort mattress and head pad
- Hand held controls:
 - Height Adjustment
 - Trendelenburg/Reverse Trendelenburg
 - Backrest Adjustment
 - Knee Flex Adjustment
 - Auto Contouring Adjustment
- Emergency Backrest Release to override electric/battery
- Battery back-up allows power functions when not plugged in

■ Standard — Not Applicable/Available

† standard rounding method used to show metric values as whole numbers



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

Distributed By:

HAYDEN MEDICAL

661-296-7200

info@haydenmedical.com

www.haydenmedical.com

www.haydenophthalmic.com

graham field

GF2300026Rev/C24