

Distributed By:

**HAYDEN MEDICAL**

661-296-7200

info@haydenmedical.com

www.haydenmedical.com

www.haydenophthalmic.com

INTRODUCING THE  
**ESC<sup>2</sup>** OPHTHALMOLOGY  
PLATFORM



DESIGNED TO BOOST  
**PRODUCTIVITY**

STREAMLINE PATIENT HANDLING  
IMPROVE PATIENT FLOW AND THROUGHPUT  
ENHANCE PATIENT, SURGEON AND STAFF SATISFACTION

**graham**  
**field**

# ESC<sup>2</sup> DESIGNED BY CLINICIANS FOR CLINICIANS

THE MOST **ADVANCED EYE CHAIR** IN TODAY'S MARKET



Patented Articulating Head Rest with Infinite Positioning

Contoured, Multi-density Pad with **DURA-GEL**® Insert for Greater Patient Comfort

Integrated Powered Leg Extension for Taller Patients

**REVERSE TRENDLENBURG** to Improve Surgical Exposure

Superior Surgeon Leg Clearance

Low Profile Base for Smaller Spaces and Greater Clearance


**SMART CASTER TECHNOLOGY** with Auto Braking and Advanced Steering for Improved Patient and Staff Safety

## ESC<sup>2</sup>

ADVANCED SINGLE PATIENT  
P L A T F O R M

**600 POUNDS**  
600 lb (272 kg)  
Weight Capacity  
EVENLY DISTRIBUTED  
Highest on the  
Market Today

# ESC<sup>2</sup> BUILT TODAY. DESIGNED FOR TOMORROW



Positioning at Your Fingertips:  
Intuitive Hand Pendant with Pre-Sets  
and 3 Memory Positions

Integrated Ergonomic  
Push Handles for  
Increased Staff Safety

Wider Patient Surface  
- 26 inch (66 cm) for  
Enhanced Patient  
Comfort and Safety



Counter Rotating Side Rails  
with Urethane Comfort Arm  
Rests for Improved Patient  
Safety and Comfort

Two-position Foot Rest For  
Patient Comfort and Safety



Industry Leading Height Range –  
20 to 40 inches (51 to 102 cm) for  
Greater Procedure Flexibility

**ESC<sup>2</sup>**  
THE NEXT GENERATION  
PLATFORM

# ESC<sup>2</sup> IMPROVED VERSATILITY, EFFICIENCY, AND COMFORT



0-15 degrees  
Seat Tilt for  
Patient Comfort  
and Safety



Optional Accessory:  
Full "U" Wrist Rest  
Item No: HSA078600



Complete the ESC<sup>2</sup>  
with the new Hausted  
Surgical Stool  
Item No: HSS500

Optional Accessory:  
O<sub>2</sub> Tank Holder  
Item No: HSA080009

### Standard Accessories:

Description	UOM
Telescoping IV Pole	1/ea
Safety Strap	1/ea
Utility Bag Hook	1/ea

### Optional Accessories:

Item No.	Description	UOM
HSA078600	Full "U" Wrist Rest	1/ea
HSA078500	Dual Wrist Rest	1/ea
HSA080003	OxyFlex II Oxygen Delivery System	1/ea

Item No.	Description	UOM
HSA008000	OxyFlex II Diffusion Trays	25/cs
H0103ESC2	Contoured Head Rests	3/set
HSA080014	Hi/Lo Foot Switch	1/ea

Item No.	Description	UOM
HSA080013	Shoulder Support Rails	1/pr
HSA080020	Foot End Push Handles	1/pr
HSA080009	O <sub>2</sub> Tank Holder	1/ea

### Specifications

#### ESC<sup>2</sup>EYESTR

	inch	cm <sup>†</sup>
Width:	32	81
Length:	83	211
Height: maximum	40	102
Height: minimum	20	51
Patient Surface Width:	26	66
Patient Surface Length: minimum	75	191
Patient Surface Length: maximum	83	211
Side Rail Length:	18.5	47
Side Rail Height:	10.25	26
Caster Diameter:	5	13
Multi-Density Sectional Pad: Dura-Gel™ insert, urethane cover	3	8
	degree <sup>††</sup>	
Trendelenburg:	15	
Reverse Trendelenburg:	6	
Seat Tilt:	0-15	
Backrest Range:	0-85	
	volts	
Dual Power:	110-240	
	lb	
Product Weight:	265	120
Weight Capacity: EVENLY DISTRIBUTED	600	272
Shipping Weight:	315	143

<sup>†</sup> standard rounding method used to show metric values to two decimal places  
<sup>††</sup> evenly distributed refer to user manual

ESC<sup>2</sup> Product Info



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.



CE Partner4U BV  
Esdoornlaan 13  
3951 DB Maarn  
The Netherlands  
www.cepartner4u.eu