

## Disposable Smoke Evacuation (DSE) Tubing

Quantity: 1 Pack of 50  
REF LA6150

**CAUTION:** U.S. Federal law restricts this device to sale by or on the order of a physician.

Latex Free  
Do Not Reuse

### Single Patient Use

### INDICATIONS FOR USE

DSE Tubing is a single-use 9" tubing designed to connect the Speculum and Smoke Evacuator Reducer.

### DIRECTIONS FOR USE

1. Attach DSE to the Smoke Evacuator Reducer
2. Slide the DSE tubing through the DSE tubing guide on your compatible instrument.

### NOTE:

Discard used tubing after each procedure. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.