



Distributed By:  
**HAYDEN MEDICAL**  
661-296-7200  
info@haydenmedical.com  
www.haydenmedical.com  
www.haydenophthalmic.com

LEADING IN  
STERILIZATION  
SINCE

**1925**

# T-Top Autoclave

## Efficiency with maximum reliability

**The only autoclave upgradeable to pre/post vacuum cycles.** A cost effective & durable closed door drying solution for medical & healthcare facilities.



Designed to Meet Your Healthcare  
Clinics' Everyday Sterilization Needs

# T-Top 10|11 Autoclave

Distributed By:  
**HAYDEN MEDICAL**  
661-296-7200  
info@haydenmedical.com  
www.haydenmedical.com  
www.haydenophthalmic.com

- ◆ Class S sterilizer (Closed door drying)
- ◆ Upgradable to pre/post vacuum (Class B)
- ◆ Clean water with every cycle! Separate water reservoirs for clean /used water
- ◆ T-Top 11 comes with automatic drain and fill kit included
- ◆ **Best price for TOP performance**
- ◆ T-Connect: Advanced remote diagnostics
- ◆ Large chamber capacity



P/N: T-Top 10/11 110V



## Innovation

### Best in class technology

- ◆ Time saver setting to set tests and sterilization cycles to start at any time
- ◆ Programming and customization functions
- ◆ Specially designed user interface to speed up and optimize daily workload



## Plug and Play

### No dedicated installation connections required

- ◆ Multicolour user-friendly touch screen with intuitive icons
- ◆ **Compact and ergonomic, designed to fit into your space**
- ◆ **Easy - just start and run programs**



## Traceability & Connectivity

### Advanced Technology that provides real-time data intelligence

- ◆ 2 USB Ports
- ◆ Wi-Fi connection & T-Connect software
- ◆ Remotely access cycle history, get notifications and alerts.
- ◆ Paperless documentation



## Reliable Performance

### Quality product that will do all you need and expect

- ◆ Easy to operate, maintain and use
- ◆ Internal memory - storing up to 500 cycles
- ◆ Advanced technology to deliver longevity and reliability

## Cloud based paperless software

T-Top devices feature the exclusive and technologically advanced T-Connect, a cloud based solution for your sterilization center.



# T-Top Sterilization Cycle Programs

	Icon	Cycle program	Load type	T-Top 10" cycle times including drying	T-Top 11" cycle times including drying
Class S		Unwrapped instruments 270°F	Unwrapped Instruments	21 min	37 min
		Wrapped Pouches 270°F	Wrapped Pouches	47 min	59 min
		Unwrapped Delicate 250°F	Unwrapped Delicate instruments	34 min	46 min
		Handpieces 270°F	Handpieces	50 min	70 min
	Icon	Cycle program	Load type	T-Top 10" cycle times including drying	T-Top 11" cycle times including drying
Class B		Unwrapped instruments 273°F	Unwrapped Instruments	24 min	38 min
		Wrapped Pouches 273°F	Wrapped Pouches	52 min	66 min
		Wrapped Pouches 250°F	Wrapped Delicate instruments	73 min	80 min
		Unwrapped Delicate 250°F	Unwrapped Delicate instruments	34 min	46 min

\* Cycle times above are an average for a standard load and could vary based on load weight, type, and external conditions

# T-Top Technical Specifications



	T-Top 10"	T-Top 11"
<b>Chamber size</b>	10 (W) x 18 (D)	11 (W) x 18 (D)
<b>Volume</b>	21L	27L
<b>Overall dimensions</b>	23" (D) x 18.1" (W) x 18" (H)	23" x 19.5" x 17.5"
<b>Max # of trays</b>	4 Trays	5 Trays
<b>Weight</b>	99 lbs.	115 lbs.
<b>Water filling</b>	Manual	Manual and Automatic

Distributed By:

**HAYDEN MEDICAL**  
 661-296-7200  
 info@haydenmedical.com  
 www.haydenmedical.com  
 www.haydenophthalmic.com

**Tuttnauer**    
 Innovation · Legacy · Partnership