



HAYDEN MEDICAL

LIPOPLASTY



Company Message



TOM OCSEK
Logistics Manager



HAYDEN CLELLAND
VP Sales & Marketing



SCOTT CLELLAND
Founder, CEO

Hayden Medical is the source for high quality German surgical instruments at reasonable prices. Why buy from the last supplier who charged your hospital, office or clinic a premium, when you can get your best deal right here with Hayden Medical.


Save money today by purchasing from Hayden Medical for all of your specialty plastic and cosmetic surgery instrument needs.

We can personalize your instruments with anything from plain text to the surgery center's logo or the surgeon's signature with our laser marking system. With over 50 years of combined surgical instrument experience, we know which instruments can fit your needs for any cosmetic or plastic surgery procedure.

We will chase down your hard-to-find items and take the guesswork out of your day. High quality German instrumentation, great service and affordable pricing is our specialty.

If you've used Hayden instruments for years or are a first time customer, you will see and feel the difference every time one of our instruments are in your hands.

RETURNED GOODS POLICY



All returns must be accompanied by a Returned Goods Authorization number, available by calling our Customer Service Department.

THE FOLLOWING INFORMATION IS NEEDED TO REQUEST AUTHORIZATION FOR A RETURN.

- REASON FOR RETURN
- DATE OF PURCHASE
- INVOICE OR PO NUMBER
- PRODUCT NUMBER AND SIZE
- QUANTITY TO RETURN

All authorized returns must be made within thirty (30) days of the invoice date and accompanied by the appropriate purchase order number and Returned Goods Authorization number.

ALLOWABLE SCHEDULE FOR AUTHORIZED RETURN.

- Products processed in error by Hayden Medical, Inc. will be credited in full
- Defective product will be credited in full
- Unless defective or shipped in error, any return accepted is subject to up to a 25% restocking fee (minimum of \$10.00)
- All returned items are subject to Hayden Medical, Inc. inspection and acceptance
- Product must be shipped prepaid

NON-RETURNABLE, OTHER THAN FOR MANUFACTURING DEFECT.

- Engraved instruments or instruments etched or labeled by the hospital, surgery center or surgeon's office
- Products damaged beyond repair or refurbishment or not in saleable condition
- Instruments returned more than thirty (30) days after invoice date
- Custom, laser marked, special or modified items
- Disposable items and Cannulas

Used products must be returned sterile wrapped or they will not be accepted.



SEND ALL RETURNS TO:
19425 Soledad Cyn. Road, #411
Santa Clarita, CA 91351



CANNULA SETS

LIPOPLASTY



Micro Aspiration Harvesting
Style 1 Micro Injector
Style 2 Micro Injector
Style 3 Micro Injector
Facial Anesthesia Needle
Snapper Lock for 20cc Syringe
Micro Transfer Hub Set

AUTOLOGOUS FAT TRANSFER SET (FACE)

ORDER ENTIRE SET
HC-FACE

PRODUCT CODE	DESCRIPTION
HC-MAN1215	Micro Aspiration Harvesting Cannula – 12ga x 15cm
HC-MIN11607	Style 1 Micro-Injector – 16ga x 7cm
HC-MIN21609	Style 2 Micro-Injector – 16ga x 9cm
HC-MIN31609	Style 3 Micro-Injector – 16ga x 9cm
HC-FAN1615	Facial Anesthesia Needle – 16ga x 15cm
HC-SL20	Snapper Lock for 20cc Syringe
NFT-SET	Micro Transfer Hub Set – ID 1.2mm, 1.4mm & 2.4mm



Tumescence Infiltration Cannula
Tri-Port Harvester
Style 1 Micro Injector
Snapper Lock for 60cc Syringe
Transfer Hub Female Luer Ends

AUTOLOGOUS FAT TRANSFER SET (BREAST)

ORDER ENTIRE SET
HC-BRST

PRODUCT CODE	DESCRIPTION
HC-BAN1215	Tumescence Infiltration Cannula – 12ga x 15cm
HC-III0315	Tri-Port Harvester – 3mm x 15cm
HC-MIN11215	Style 1 Micro-Injector – 12ga x 15cm
HC-SL60	Snapper Lock for 60cc Syringe
NFT-2.4	Transfer Hub, Female Luer Ends – ID 2.4mm



Mercedes Cannula (2)
Style 1 Injector (2)
Tumescence Infiltration Cannula

AUTOLOGOUS FAT TRANSFER SET (GLUTEAL)

ORDER ENTIRE SET
HC-GLUT

PRODUCT CODE	DESCRIPTION
HC-MER0420	Mercedes Cannula – 4mm x 20cm
HC-MER0425	Mercedes Cannula – 4mm x 25cm
HC-MIN10315	Style 1 Injector – 3mm x 15cm
HC-MIN10320	Style 1 Injector – 3mm x 20cm
HC-BAN1225	Tumescence Infiltration Cannula – 12ga x 25cm



Tumescence Infiltration Cannula (2)

Black Infusion Handle

TUMESCENT INFILTRATION SET

ORDER ENTIRE SET
HC-TINF

PRODUCT CODE	DESCRIPTION
HC-BAN1232	Tumescence Infiltration Cannula – 12ga x 32cm
HC-BAN1415	Tumescence Infiltration Cannula – 14ga x 15cm
HC-INFU-C-BLK	Infusion Handle with On/Off Control – Black



CANNULA HO

LIPOPLASTY

Hayden Medical's cannulas and cannula accessories are made of high-quality, medical-grade stainless steel and made to order in the USA. If there is a size or pattern that you don't see, we can make it. Please contact us directly with custom cannula requests.



HAYDEN MULTI-HOLE HARVESTER

270° of extraction and 18 holes for low-volume liposuction procedures (available with angled shaft). Typically used to harvest for small-volume fat grafting to areas like the face and hands.

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-SAT0315	3mm	15cm	LL / LH / TM
HC-SAT0326	3mm	26cm	LL / LH / TM
HC-SAT0332	3mm	32cm	LL / LH / TM

MERCEDES

360° of extraction and 3 holes for high-volume liposuction procedures. Typically the most popular and commonly used liposuction cannula.

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-MER0326	3mm	26cm	LL / LH / TM
HC-MER0332	3mm	32cm	LL / LH / TM
HC-MER0426	4mm	26cm	LL / LH / TM
HC-MER0432	4mm	32cm	LL / LH / TM
HC-MER0526	5mm	26cm	LL / LH / TM
HC-MER0532	5mm	32cm	LL / LH / TM

BASKET

360° of extraction and 3 flared openings for large-volume liposuction procedures. These can be referred to as a "Flared Mercedes" cannula as well.

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-BAS0326	3mm	26cm	LL / LH / TM
HC-BAS0332	3mm	32cm	LL / LH / TM
HC-BAS0426	4mm	26cm	LL / LH / TM
HC-BAS0432	4mm	32cm	LL / LH / TM
HC-BAS0526	5mm	26cm	LL / LH / TM
HC-BAS0532	5mm	32cm	LL / LH / TM

MERCEDES PLUS (DOUBLE MERCEDES)

360° of extraction and 6 holes for high volume liposuction procedures.

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-MER+0326	3mm	26cm	LL / LH / TM
HC-MER+0332	3mm	32cm	LL / LH / TM
HC-MER+0426	4mm	26cm	LL / LH / TM
HC-MER+0432	4mm	32cm	LL / LH / TM
HC-MER+0526	5mm	26cm	LL / LH / TM
HC-MER+0532	5mm	32cm	LL / LH / TM

CANDY CANE

360° of extraction for medium to high-volume liposuction procedures.

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-CC0326	3mm	26cm	LL / LH / TM
HC-CC0332	3mm	32cm	LL / LH / TM
HC-CC0426	4mm	26cm	LL / LH / TM
HC-CC0432	4mm	32cm	LL / LH / TM
HC-CC0526	5mm	26cm	LL / LH / TM
HC-CC0532	5mm	32cm	LL / LH / TM

LE PATTERNS

LIPOPLASTY



SPATULA

90° of extraction and a single hole for low to medium-volume liposuction procedures. Non-traumatic spatula tips are ideal for delicate areas like the face and neck.

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-SPA0215	2mm	15cm	LL / LH / TM
HC-SPA0315	3mm	15cm	LL / LH / TM
HC-SPA0415	4mm	15cm	LL / LH / TM
HC-SPA0515	5mm	15cm	LL / LH / TM

HUB OPTIONS



SPATULA 2-PORT

90° of extraction and 2 holes for low to medium-volume liposuction procedures. Non-traumatic spatula tips are ideal for delicate areas like the face and neck.

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-SP20215	2mm	15cm	LL / LH / TM
HC-SP20315	3mm	15cm	LL / LH / TM
HC-SP20415	4mm	15cm	LL / LH / TM
HC-SP20515	5mm	15cm	LL / LH / TM



MODERN SPIRAL

270° of extraction with multiple holes for high-volume aspiration. The Modern Spiral cannula is becoming the new favorite of surgeons who traditionally use a Tri-Port or Mercedes cannula.

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-MS0326	3mm	26cm	LL / LH / TM
HC-MS0332	3mm	32cm	LL / LH / TM
HC-MS0426	4mm	26cm	LL / LH / TM
HC-MS0432	4mm	32cm	LL / LH / TM
HC-MS0526	5mm	26cm	LL / LH / TM
HC-MS0532	5mm	32cm	LL / LH / TM



TRI-PORT

270° of extraction and 3 holes for high-volume liposuction procedures

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-III0215	2mm	15cm	LL / LH / TM
HC-III0226	2mm	26cm	LL / LH / TM
HC-III0315	3mm	15cm	LL / LH / TM
HC-III0326	3mm	26cm	LL / LH / TM
HC-III0332	3mm	32cm	LL / LH / TM
HC-III0415	4mm	15cm	LL / LH / TM
HC-III0426	4mm	26cm	LL / LH / TM
HC-III0432	4mm	32cm	LL / LH / TM
HC-III0526	5mm	26cm	LL / LH / TM
HC-III0532	5mm	32cm	LL / LH / TM



TRI-PORT PLUS FOR LARGE-VOLUME HARVESTING

270° of extraction and 6 holes for high-volume liposuction procedures.

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-III+0326	3mm	26cm	LL / LH / TM
HC-III+0332	3mm	32cm	LL / LH / TM
HC-III+0426	4mm	26cm	LL / LH / TM
HC-III+0432	4mm	32cm	LL / LH / TM
HC-III+0526	5mm	26cm	LL / LH / TM
HC-III+0532	5mm	32cm	LL / LH / TM



INJECTION & INFILT

MICRO INJECTION



STYLE 1

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-MIN11607	16ga	7cm	LUER LOCK
HC-MIN11609	16ga	9cm	LUER LOCK
HC-MIN11807	18ga	7cm	LUER LOCK
HC-MIN11809	18ga	9cm	LUER LOCK
HC-MIN11905	19ga	5cm	LUER LOCK
HC-MIN12003	20ga	3cm	LUER LOCK
HC-MIN12005	20ga	5cm	LUER LOCK



STYLE 2

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-MIN21607	16ga	7cm	LUER LOCK
HC-MIN21609	16ga	9cm	LUER LOCK
HC-MIN21807	18ga	7cm	LUER LOCK
HC-MIN21809	18ga	9cm	LUER LOCK



STYLE 3

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-MIN31607	16ga	7cm	LUER LOCK
HC-MIN31609	16ga	9cm	LUER LOCK
HC-MIN31807	18ga	7cm	LUER LOCK
HC-MIN31809	18ga	9cm	LUER LOCK



"V" DISSECTOR

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-MINV1215	12ga	15cm	LUER LOCK
HC-MINV1225	12ga	25cm	LUER LOCK
HC-MINV1415	14ga	15cm	LUER LOCK
HC-MINV1409	14ga	9cm	LUER LOCK
HC-MINV1607	16ga	7cm	LUER LOCK
HC-MINV1609	16ga	9cm	LUER LOCK
HC-MINV1807	18ga	7cm	LUER LOCK

MICRO ASPIRATION



MERCEDES

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-MER1215	12ga	15cm	LUER LOCK
HC-MER1225	12ga	25cm	LUER LOCK
HC-MER1415	14ga	15cm	LUER LOCK
HC-MER1425	14ga	25cm	LUER LOCK
HC-MER0315	3mm	15cm	LUER LOCK
HC-MER0325	3mm	25cm	LUER LOCK



FAT ASPIRATION CANNULA (FAC)

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-MAN1215	12ga	15cm	LUER LOCK
HC-MAN1220	12ga	20cm	LUER LOCK
HC-MAN1225	12ga	25cm	LUER LOCK
HC-MAN1415	14ga	15cm	LUER LOCK
HC-MAN1425	14ga	25cm	LUER LOCK

RATION CANNULAS

BODY TISSUE INJECTION

TUMESCENT INFILTRATION



STYLE 1

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-MIN11215	12ga	15cm	LUER LOCK
HC-MIN11220	12ga	20cm	LUER LOCK
HC-MIN11225	12ga	25cm	LUER LOCK
HC-MIN11415	14ga	15cm	LUER LOCK
HC-MIN11420	14ga	20cm	LUER LOCK
HC-MIN11425	14ga	25cm	LUER LOCK




BODY ANESTHESIA NEEDLE

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-BAN1015	10ga	15cm	LUER LOCK
HC-BAN1032	10ga	32cm	LUER LOCK
HC-BAN1215	12ga	15cm	LUER LOCK
HC-BAN1225	12ga	25cm	LUER LOCK
HC-BAN1232	12ga	32cm	LUER LOCK
HC-BAN1415	14ga	15cm	LUER LOCK
HC-BAN1425	14ga	25cm	LUER LOCK
HC-BAN1432	14ga	32cm	LUER LOCK
HC-BAN1615	16ga	15cm	LUER LOCK

TISSUE TRANSFER

SOFT TISSUE TRANSFER SYSTEM

SYRINGE, LUER LOCK TO LUER LOCK




PRODUCT CODE	DIAMETER
HC-LL3	3mm
HC-LL4	4mm




FACIAL ANESTHESIA NEEDLE

PRODUCT CODE	DIAMETER	LENGTH	HUB OPTIONS
HC-FAN1415	14ga	15cm	LUER LOCK
HC-FAN1615	16ga	15cm	LUER LOCK

MICRO TRANSFER HUB SET



1.2mm



1.4mm



2.4mm

PRODUCT CODE	DESCRIPTION
NFT-SET	Set of 3



TWIST LOCK

PRODUCT CODE	DESCRIPTION
HC-MJ-TL60	for 60cc Monoject Syringe
HC-BD-TL60	for 60cc BD Syringe
HC-MJ-TL35	for 35cc Monoject Syringe
HC-BT-TL20	for 20cc BD and 20cc Terumo Syringe
HC-BD-TL30	for 30cc BD Syringe
HC-BD-TL10	for 10cc BD Syringe
HC-MJ-TL10	for 10cc Monoject and 12cc Monoject Syringe



LIPOSUCTION HANDLES & ACCESSORIES

INFUSION HANDLE 1" with On/Off Control



PRODUCT CODE	DESCRIPTION
HC-INFU-C-BLK	A flow on/off control infusion handle for use with roller pump devices.


INFUSION HANDLE



PRODUCT CODE	DESCRIPTION
HC-INFU-NC-L	Standard Infusion handle, No Control, Ideal for infusion needle 1" diameter.

PLASTIC STERILE TRAYS

- Autoclavable
- Silicone Finger-Tip Mats



PRODUCT CODE	DESCRIPTION
030-105	4" x 7-1/2" x 3/4"
030-105.5	4" x 7-1/2" x 1-1/2"
030-107	6" x 10" x 3/4"
030-108	6" x 10" x 1-1/2"
030-108D	6" x 10" x 1-1/2" Double-Decker
030-111	10" x 15" x 3/4"
030-112	10" x 15" x 1-1/2"

NEEDLE CONNECTOR



PRODUCT CODE	DESCRIPTION
HC-NC	Needle connector for use with tubing

ASPIRATION HANDLES



PRODUCT CODE	DESCRIPTION
HC-AH-S	Standard Aspiration Handle 3/4"
HC-AH-L	Standard Aspiration Handle 1"

AWLS



PRODUCT CODE	DESCRIPTION
HC-EP3	3mm
HC-EP4	4mm
HC-EP5	5mm

SNAPPER LOCKS



PRODUCT CODE	DESCRIPTION
HC-BT-SL10	for 10cc BD and 10cc Terumo Syringe
HC-MJ-SL10	for 10cc Monoject and 12cc Monoject Syringe
HC-BD-SL20	for 20cc BD Syringe
HC-MJ-SL20	for 20cc Monoject Syringe
HC-BD-SL30	for 30cc BD Syringe
HC-MJ-SL35	for 35cc Monoject Syringe
HC-TR-SL60	for 60cc Terumo Syringe
HC-MB-SL60	for 60cc Monoject and B-D Syringe

ADIPOSE TISSUE PROCESSING

PRISM PARTICLE SIZER™

The **PRISM PARTICLE SIZER™** is a particle-sizing and filtering device for human tissue and/or fluid. When used to process fat aspirate, the device creates precision-sized micro-fragmented adipose tissue (MFAT) across a selectable dynamic sizing range of 3000-µm to 200-µm.

BENEFITS OF PPS DEVICE INCLUDES:

- More ergonomic single pass sizing & filtering
- Creates the most purified processed fat aspirate
- Cost-effective, reusable, Class I medical device
- Patent pending and manufactured in the USA
- Recommended PPS device accessories include:

PRISM CAROUSEL™

PRISM TRANSFER CANNULA™

PRISM CANNULA CLEANER™

PRODUCT CODE	DESCRIPTION					
PPS-Face	2000-µm	1500-µm	900-µm	500-µm	300-µm	
PPS-Breast	2500-µm	1500-µm	1200-µm			
PPS-Hands	2000-µm	1200-µm	1000-µm			
PPS-Buttocks	3000-µm	2500-µm	2000-µm	1500-µm		
PPS-Regenerative	3000-µm	2500-µm	2000-µm	1500-µm		
PPS-Set	3000-µm	2500-µm	2000-µm	1500-µm	1200-µm	
	1000-µm	900-µm	800-µm	700-µm	600-µm	
	500-µm	400-µm	300-µm	200-µm		



PRISM PARTICLE SIZER ACCESSORIES

PRISM CAROUSEL™




PRODUCT CODE	DESCRIPTION
P-CAR-6	One Decanting/Storage Rack for 6 Luer-Lock Syringes
P-CAR-12	One Decanting/Storage Rack for 12 Luer-Lock Syringes

PRISM CANNULA CLEANER™



PRODUCT CODE	DESCRIPTION
PCC	One reusable nylon tool for cleaning and removing debris when reprocessing PPS device

PRISM TRANSFER CANNULA™



PRODUCT CODE	DESCRIPTION
PTC-SET	Set of one shorter and one longer Cannula for transferring fat aspirate while using the PPS device system

The **PRISM PARTICLE SIZER™** medical device and its recommended accessories are FDA-registered and manufactured by HOGUE SURGICAL (Maple Grove, MN, USA)

PUMPS & TUBING

INFILTRATION PUMP

PRODUCT CODE

HC-IP

- Delivers High-Volume of Tumescant Fluid
- Can Sit on Aspirator
- Compact and Minimal Noise



ASPIRATOR PUMP

PRODUCT CODE

HC-AP

- High-Torque Suction
- Portable and Compact
- Integrated Pressure Gauge
- Minimal Device Noise



INFILTRATION TUBING



LATEX FREE
SINGLE USE STERILE

PRODUCT CODE

HC-IDT

DESCRIPTION

Infiltration pump tubing
double spike tubing, box 10



LATEX FREE
SINGLE USE STERILE

PRODUCT CODE

HC-IST

DESCRIPTION

Infiltration pump tubing, box 10

★ **MOST POPULAR**

ASPIRATION TUBING



LATEX FREE
SINGLE USE STERILE

PRODUCT CODE

HC-AT

DESCRIPTION

Flexible ribbed aspiration tubing,
1/4" ID, 12 ft length, box 20

ASPIRATION CANISTER

*Please note canister and lid/liner sold separatley



HC-L10

HC-C01

PRODUCT CODE

HC-C01

HC-L10

DESCRIPTION

Canister, holds up to 1500cc

Liner, soft-shell plastic with snap
top lid holds up to 1500cc, box 10

LASER & FACIAL PLASTIC SURGERY

OCULAR LASER SHIELD WITH SUCTION CUP



PRODUCT CODE	DESCRIPTION			
OP21-261P	X-Small	22.0 x 19.5mm	Pair	
OP21-262P	Small	24.0 x 21.0mm	Pair	
OP21-264P	Medium	25.5 x 22.0mm	Pair	
OP21-266P	Large	26.5 x 23.0mm	Pair	

PLASTIC OCULAR SHIELD WITH SUCTION CUP



PRODUCT CODE	DESCRIPTION			
OP23-201BP	X-Small	23.5 x 21.5mm	Pair	
OP23-202BP	Small	26.0 x 23.5mm	Pair	
OP23-204BP	Medium	27.5 x 24.5mm	Pair	
OP23-206BP	Large	28.5 x 25.5mm	Pair	

OCULAR LASER SHIELD WITH OFFSET HANDLE



PRODUCT CODE	DESCRIPTION			
OP21-261HP	X-Small	22.0 x 19.5mm	Pair	
OP21-262HP	Small	24.0 x 21.0mm	Pair	
OP21-264HP	Medium	25.5 x 22.0mm	Pair	
OP21-266HP	Large	26.5 x 23.0mm	Pair	

PLASTIC OCULAR SHIELD WITH HANDLE



PRODUCT CODE	DESCRIPTION			
OP23-201BHP	X-Small	23.5 x 21.5mm	Pair	
OP23-202BHP	Small	26.0 x 23.5mm	Pair	
OP23-204BHP	Medium	27.5 x 24.5mm	Pair	
OP23-206BHP	Large	28.5 x 25.5mm	Pair	


DESMARRES RETRACTOR NON-CONDUCTIVE

210-198NC



PRODUCT CODE	DESCRIPTION
210-194NC	11mm wide
210-195NC	13mm wide
210-196NC	15mm wide
210-197NC	17mm wide
210-198NC	8 / 10mm wide - double ended

CASTROVIEJO NEEDLE HOLDER 5-1/2" TUNGSTEN



PRODUCT CODE	DESCRIPTION
091-828TC	Straight, w/ Lock, Spring Handle - Smooth
091-829TC	Straight, w/ Lock, Spring Handle - Serrated
091-832TC	Curved, w/ Lock, Spring Handle - Smooth
091-833TC	Curved, w/ Lock, Spring Handle - Serrated

HAYDEN BLEPHAROPLASTY SCISSORS 4-1/2" TUNGSTEN



PRODUCT CODE	DESCRIPTION
500-960	Curved, flat tipped serrated blade

CASTROVIEJO SUTURING FORCEPS 4"



PRODUCT CODE	DESCRIPTION
120-952D	0.5mm, Diamond Dust Platform
120-953	0.12mm, w/ Tying Platform
120-954	0.3mm, w/ Tying Platform
120-955	0.5mm, w/ Tying Platform



SURGICAL SUITE EQUIPMENT

Outfit your surgical suite and procedure rooms with our wide range of new and refurbished equipment. We offer major brands and can match your needs with multiple choices to suit your specific budget. We pride ourselves on the quality of our products, excellent customer service and service after the sale.

• SURGICAL LIGHTING

QALED070SC



QALED0100DC



• STAINLESS STEEL TABLES AND STANDS

HCM731



HCM513



HCM520



• PROCEDURE CHAIRS

ESCEYE-ST500



• PATIENT MONITORS



• TABLE TOP AUTOCLAVES



• DOCTOR'S STOOL

EDS-1060G



• SURGICAL STOOL

S900-W



ANESTHESIA MACHINES, AUTOCLAVES, SURGICAL C-ARMS, C-ARM TABLES, DEFIBRILLATORS, ELECTROSURGICAL UNITS, ENDOSCOPY EQUIPMENT, EXAM TABLES, MEDICAL STAINLESS STEEL, MEDICAL STRETCHERS, PATIENT MONITORS, PROCEDURE TABLES, SURGICAL MICROSCOPES, SURGICAL LIGHTING, SURGICAL TABLES, WARMING CABINETS

FACE, BREAST & BODY

Contact us to discuss our extensive line of instruments for cosmetic face, breast and body procedures. The breast instruments below help increase safety and minimize surgical time.

• LIGHT SOURCE, LIGHT CABLE & BREAST RETRACTOR



• MONOPOLAR CAUTERY FORCEPS



BIPOLAR FORCEPS ALSO AVAILABLE

• ELECTROSURGICAL UNIT & SMOKE EVACUATOR

