




STERILIZATION CONTAINERS



**HAYDEN
MEDICAL^{INC}**

(877) 296-7205 tel
(877) 296-7206 fax



	1/1 Containers and accessories	2-5
	3/4 Containers and accessories	6-9
	1/2 Containers and accessories	10-13
	Flat half-size-containers and accessories	14
	US Containers and accessories	15
	Small containers and accessories	16
	Dental Container and accessories	17
	Optic Container	18
	Standard accessories	19-21
	Sheet baskets and covers	22
	Sterilizing cases	23
	Wire baskets / Container for endoscopes	24
	Storage- and sterilization rack / silicone bar	25

Container models and dimensions	Safety lid		Safety lid	
	Perforated lid Unperforated bottom	Perforated lid Perforated bottom	Perforated lid Unperforated bottom	Perforated lid Perforated bottom
580 x 280 x 100 mm 	R 110.10 R 120.10 R 130.10 R 140.10 R 150.10 R 160.10	R 111.10 R 121.10 R 131.10 R 141.10 R 151.10 R 161.10	R 112.10 R 122.10 R 132.10 R 142.10 R 152.10 R 162.10	R 113.10 R 123.10 R 133.10 R 143.10 R 153.10 R 163.10
580 x 280 x 135 mm 	R 110.13 R 120.13 R 130.13 R 140.13 R 150.13 R 160.13	R 111.13 R 121.13 R 131.13 R 141.13 R 151.13 R 161.13	R 112.13 R 122.13 R 132.13 R 142.13 R 152.13 R 162.13	R 113.13 R 123.13 R 133.13 R 143.13 R 153.13 R 163.13
580 x 280 x 150 mm 	R 110.15 R 120.15 R 130.15 R 140.15 R 150.15 R 160.15	R 111.15 R 121.15 R 131.15 R 141.15 R 151.15 R 161.15	R 112.15 R 122.15 R 132.15 R 142.15 R 152.15 R 162.15	R 113.15 R 123.15 R 133.15 R 143.15 R 153.15 R 163.15
580 x 280 x 200 mm 	R 110.20 R 120.20 R 130.20 R 140.20 R 150.20 R 160.20	R 111.20 R 121.20 R 131.20 R 141.20 R 151.20 R 161.20	R 112.20 R 122.20 R 132.20 R 142.20 R 152.20 R 162.20	R 113.20 R 123.20 R 133.20 R 143.20 R 153.20 R 163.20
580 x 280 x 260 mm 	R 110.26 R 120.26 R 130.26 R 140.26 R 150.26 R 160.26	R 111.26 R 121.26 R 131.26 R 141.26 R 151.26 R 161.26	R 112.26 R 122.26 R 132.26 R 142.26 R 152.26 R 162.26	R 113.26 R 123.26 R 133.26 R 143.26 R 153.26 R 163.26

Unperforated bottom



580 x 280 x 100 mm **R 105.10**
 580 x 280 x 135 mm **R 105.13**
 580 x 280 x 150 mm **R 105.15**
 580 x 280 x 200 mm **R 105.20**
 580 x 280 x 260 mm **R 105.26**

Perforated bottom



580 x 280 x 100 mm **R 106.10**
 580 x 280 x 135 mm **R 106.13**
 580 x 280 x 150 mm **R 106.15**
 580 x 280 x 200 mm **R 106.20**
 580 x 280 x 260 mm **R 106.26**

580 x 280 mm

Perforated lid

**R 100.01****R 100.02****R 100.03****R 100.04****R 100.05****R 100.06**

580 x 280 mm

Safety lid

**R 101.01****R 101.02****R 101.03****R 101.04****R 101.05****R 101.06**

Wire baskets



485 x 250 x 30 mm **R 190.03**
 485 x 250 x 50 mm **R 190.05**
 485 x 250 x 70 mm **R 190.07**
 485 x 250 x 100 mm **R 190.10**

Wire baskets



540 x 250 x 30 mm **R 191.03**
 540 x 250 x 50 mm **R 191.05**
 540 x 250 x 70 mm **R 191.07**
 540 x 250 x 100 mm **R 191.10**

Wire base



545 x 255 x 30 mm **R 193.00**

Silicon mat



470 x 230 mm **R 852.40**
 520 x 230 mm **R 852.30**

Tension clips for towels



545 x 255 mm **R 195.00**

Drape




130 x 60 cm **R 095.07**
 140 x 100 cm **R 095.08**



Perforated lid
Unperforated bottom



580 x 280 x 115 mm



R 110.10 BB






R 120.10 BB

R 130.10 BB


R 140.10 BB

R 150.10 BB

R 160.10 BB



580 x 280 x 140 mm



R 110.13 BB

R 120.13 BB


R 130.13 BB

R 140.13 BB

R 150.13 BB

R 160.13 BB

580 x 280 x 155 mm



R 110.15 BB

R 120.15 BB


R 130.15 BB

R 140.15 BB

R 150.15 BB

R 160.15 BB

580 x 280 x 205 mm



R 110.20 BB

R 120.20 BB


R 130.20 BB

R 140.20 BB

R 150.20 BB

R 160.20 BB

580 x 280 x 265 mm



R 110.26 BB

R 120.26 BB

R 130.26 BB

R 140.26 BB

R 150.26 BB

R 160.26 BB


1. Before sterilization, the valves are in the relaxed position.
2. During sterilization, the vacuum completely opens the valves and the dry air is released from the container.
3. At the end of the sterilization cycle, the vacuum completely opens the valves in the opposite direction.
4. At the end of the sterilization cycle, when the pressure level is equalized, the valves are closed.
5. After sterilization, the valves are in the relaxed position.



Perforated lid
Perforated bottom



580 x 280 x 135 mm



R 111.10 BB

R 121.10 BB

R 131.10 BB

R 141.10 BB

R 151.10 BB

R 161.10 BB






1. Before sterilization, the valves are in the relaxed position.

2. During sterilization, the vacuum completely opens the valves and the dry air is released from the container.


3. At the end of the sterilization cycle, the vacuum completely opens the valves in the opposite direction.

4. At the end of the sterilization cycle, when the pressure level is equalized, the valves are closed.

5. After sterilization, the valves are in the relaxed position.



580 x 280 x 160 mm



R 111.13 BB

R 121.13 BB


R 131.13 BB

R 141.13 BB

R 151.13 BB

R 161.13 BB

580 x 280 x 175 mm



R 111.15 BB

R 121.15 BB


R 131.15 BB

R 141.15 BB

R 151.15 BB

R 161.15 BB

580 x 280 x 235 mm



R 111.20 BB

R 121.20 BB


R 131.20 BB

R 141.20 BB

R 151.20 BB

R 161.20 BB

580 x 280 x 295 mm



R 111.26 BB

R 121.26 BB

R 131.26 BB

R 141.26 BB

R 151.26 BB

R 161.26 BB

Container models and dimensions	Safety lid		Safety lid	
	Perforated lid Unperforated bottom	Perforated lid Perforated bottom	Perforated lid Unperforated bottom	Perforated lid Perforated bottom
465 x 280 x 100 mm 	R 210.10 R 220.10 R 230.10 R 240.10 R 250.10 R 260.10	R 211.10 R 221.10 R 231.10 R 241.10 R 251.10 R 261.10	R 212.10 R 222.10 R 232.10 R 242.10 R 252.10 R 262.10	R 213.10 R 223.10 R 233.10 R 243.10 R 253.10 R 263.10
465 x 280 x 135 mm 	R 210.13 R 220.13 R 230.13 R 240.13 R 250.13 R 260.13	R 211.13 R 221.13 R 231.13 R 241.13 R 251.13 R 261.13	R 212.13 R 222.13 R 232.13 R 242.13 R 252.13 R 262.13	R 213.13 R 223.13 R 233.13 R 243.13 R 253.13 R 263.13
465 x 280 x 150 mm 	R 210.15 R 220.15 R 230.15 R 240.15 R 250.15 R 260.15	R 211.15 R 221.15 R 231.15 R 241.15 R 251.15 R 261.15	R 212.15 R 222.15 R 232.15 R 242.15 R 252.15 R 262.15	R 213.15 R 223.15 R 233.15 R 243.15 R 253.15 R 263.15

Unperforated bottom



465 x 280 x 100 mm **R 205.10**
465 x 280 x 135 mm **R 205.13**
465 x 280 x 150 mm **R 205.15**

Perforated bottom




465 x 280 x 100 mm **R 206.10**
465 x 280 x 135 mm **R 206.13**
465 x 280 x 150 mm **R 206.15**

465 x 280 mm




Perforated lid



465 x 280 mm



Safety lid



Wire baskets



405 x 250 x 30 mm **R 290.03**
405 x 250 x 50 mm **R 290.05**
405 x 250 x 70 mm **R 290.07**
405 x 250 x 100 mm **R 290.10**

Wire base



405 x 250 x 30 mm **R 293.03**

Silicon mat



380 x 230 mm **R 852.10**
440 x 230 mm **R 852.20**

Tension clips for towels












405 x 255 mm **R 295.00**

Drape



140 x 100 cm **R 095.08**


	Perforated lid Unperforated bottom 	
465 x 280 x 115 mm 	R 210.10 BB R 220.10 BB R 230.10 BB R 240.10 BB R 250.10 BB R 260.10 BB	- Before sterilization, the valves are in the relaxed position. - During sterilization, the vacuum completely opens the valves and the dry air is released from the container.  
465 x 280 x 140 mm 	R 210.13 BB R 220.13 BB R 230.13 BB R 240.13 BB R 250.13 BB R 260.13 BB	- At the end of the sterilization cycle, the vacuum completely opens the valves in the opposite direction. - At the end of the sterilization cycle, when the pressure level is equalized, the valves are closed. - After sterilization, the valves are in the relaxed position.  
465 x 280 x 155 mm 	R 210.15 BB R 220.15 BB R 230.15 BB R 240.15 BB R 250.15 BB R 260.15 BB	



Perforated lid
Perforated bottom



465 x 280 x 135 mm



R 211.10 BB

R 221.10 BB

R 231.10 BB

R 241.10 BB

R 251.10 BB

R 261.10 BB






1. Before sterilization, the valves are in the relaxed position.

2. During sterilization, the vacuum completely opens the valves and the dry air is released from the container.


3. At the end of the sterilization cycle, the vacuum completely opens the valves in the opposite direction.

4. At the end of the sterilization cycle, when the pressure level is equalized, the valves are closed.

5. After sterilization, the valves are in the relaxed position.



465 x 280 x 160 mm



R 211.13 BB

R 221.13 BB


R 231.13 BB

R 241.13 BB

R 251.13 BB

R 261.13 BB

465 x 280 x 175 mm



R 211.15 BB


R 221.15 BB

R 231.15 BB

R 241.15 BB

R 251.15 BB

R 261.15 BB

Container models and dimensions	Safety lid		Safety lid	
	Perforated lid Unperforated bottom	Perforated lid Perforated bottom	Perforated lid Unperforated bottom	Perforated lid Perforated bottom
285 x 280 x 100 mm 	R 310.10 R 320.10 R 330.10 R 340.10 R 350.10 R 360.10	R 311.10 R 321.10 R 331.10 R 341.10 R 351.10 R 361.10	R 312.10 R 322.10 R 332.10 R 342.10 R 352.10 R 362.10	R 313.10 R 323.10 R 333.10 R 343.10 R 353.10 R 363.10
285 x 280 x 135 mm 	R 310.13 R 320.13 R 330.13 R 340.13 R 350.13 R 360.13	R 311.13 R 321.13 R 331.13 R 341.13 R 351.13 R 361.13	R 312.13 R 322.13 R 332.13 R 342.13 R 352.13 R 362.13	R 313.13 R 323.13 R 333.13 R 343.13 R 353.13 R 363.13
285 x 280 x 150 mm 	R 310.15 R 320.15 R 330.15 R 340.15 R 350.15 R 360.15	R 311.15 R 321.15 R 331.15 R 341.15 R 351.15 R 361.15	R 312.15 R 322.15 R 332.15 R 342.15 R 352.15 R 362.15	R 313.15 R 323.15 R 333.15 R 343.15 R 353.15 R 363.15
285 x 280 x 200 mm 	R 310.20 R 320.20 R 330.20 R 340.20 R 350.20 R 360.20	R 311.20 R 321.20 R 331.20 R 341.20 R 351.20 R 361.20	R 312.20 R 322.20 R 332.20 R 342.20 R 352.20 R 362.20	R 313.20 R 323.20 R 333.20 R 343.20 R 353.20 R 363.20
285 x 280 x 260 mm 	R 310.26 R 320.26 R 330.26 R 340.26 R 350.26 R 360.26	R 311.26 R 321.26 R 331.26 R 341.26 R 351.26 R 361.26	R 312.26 R 322.26 R 332.26 R 342.26 R 352.26 R 362.26	R 313.26 R 323.26 R 333.26 R 343.26 R 353.26 R 363.26

Unperforated bottom



280 x 280 x 100 mm **R 305.10**
 280 x 280 x 135 mm **R 305.13**
 280 x 280 x 150 mm **R 305.15**
 280 x 280 x 200 mm **R 305.20**
 280 x 280 x 260 mm **R 305.26**

Perforated bottom



280 x 280 x 100 mm **R 306.10**
 280 x 280 x 135 mm **R 306.13**
 280 x 280 x 150 mm **R 306.15**
 280 x 280 x 200 mm **R 306.20**
 280 x 280 x 260 mm **R 306.26**

285 x 280 mm

Perforated lid



R 300.01
R 300.02
R 300.03
R 300.04
R 300.05
R 300.06

285 x 280 mm

Safety lid



R 301.01
R 301.02
R 301.03
R 301.04
R 301.05
R 301.06

Wire baskets



255 x 250 x 30 mm **R 390.03**
 255 x 250 x 50 mm **R 390.05**
 255 x 250 x 70 mm **R 390.07**
 255 x 250 x 100 mm **R 390.10**

Wire base



255 x 255 x 30 mm **R 393.00**

Silicon mat



250 x 235 mm **R 851.00**

Tension clips for towels



255 x 255 mm **R 395.00**


Drape




90 x 60 cm **R 095.04**
 90 x 80 cm **R 095.05**
 130 x 45 cm **R 095.06**
 130 x 60 cm **R 095.07**



Perforated lid
Unperforated bottom



285 x 280 x 115 mm



R 310.10 BB

R 320.10 BB

R 330.10 BB

R 340.10 BB

R 350.10 BB

R 360.10 BB






1. Before sterilization, the valves are in the relaxed position.

2. During sterilization, the vacuum completely opens the valves and the dry air is released from the container.


3. At the end of the sterilization cycle, the vacuum completely opens the valves in the opposite direction.

4. At the end of the sterilization cycle, when the pressure level is equalized, the valves are closed.

5. After sterilization, the valves are in the relaxed position.



285 x 280 x 140 mm



R 310.13 BB

R 320.13 BB


R 330.13 BB

R 340.13 BB

R 350.13 BB

R 360.13 BB

285 x 280 x 155 mm



R 310.15 BB

R 320.15 BB


R 330.15 BB

R 340.15 BB

R 350.15 BB

R 360.15 BB

285 x 280 x 205 mm



R 310.20 BB

R 320.20 BB


R 330.20 BB

R 340.20 BB

R 350.20 BB

R 360.20 BB

285 x 280 x 265 mm



R 310.26 BB

R 320.26 BB

R 330.26 BB

R 340.26 BB

R 350.26 BB


R 360.26 BB



Perforated lid
Perforated bottom



285 x 280 x 135 mm



R 311.10 BB

R 321.10 BB

R 331.10 BB

R 341.10 BB

R 351.10 BB

R 361.10 BB






1. Before sterilization, the valves are in the relaxed position.

2. During sterilization, the vacuum completely opens the valves and the dry air is released from the container.


3. At the end of the sterilization cycle, the vacuum completely opens the valves in the opposite direction.

4. At the end of the sterilization cycle, when the pressure level is equalized, the valves are closed.

5. After sterilization, the valves are in the relaxed position.



285 x 280 x 160 mm



R 311.13 BB

R 321.13 BB


R 331.13 BB

R 341.13 BB

R 351.13 BB

R 361.13 BB

285 x 280 x 175 mm



R 311.15 BB

R 321.15 BB


R 331.15 BB

R 341.15 BB

R 351.15 BB

R 361.15 BB

285 x 280 x 235 mm



R 311.20 BB

R 321.20 BB


R 331.20 BB

R 341.20 BB

R 351.20 BB

R 361.20 BB

285 x 280 x 295 mm



R 311.26 BB

R 321.26 BB

R 331.26 BB

R 341.26 BB

R 351.26 BB

R 361.26 BB

Container models and dimensions

	Perforated lid Unperforated bottom	Perforated lid Perforated bottom	Unperforated bottom	Perforated bottom
285 x 280 x 55 mm 	R 310.05 R 320.05 R 330.05 R 340.05 R 350.05 R 360.05		R 305.05 	
285 x 280 x 85 mm 	R 310.08 R 320.08 R 330.08 R 340.08 R 350.08 R 360.08	R 311.08 R 321.08 R 331.08 R 341.08 R 351.08 R 361.08	R 305.08 	R 306.08 
285 x 280 mm 	R 304.01 R 304.02 R 304.03 R 304.04 R 304.05 R 304.06	Perforated lid	Drape 	90 x 60 cm R 095.04 90 x 80 cm R 095.05
Wire baskets 	255 x 250 x 30 mm R 390.03 255 x 250 x 50 mm R 390.05	Wire base 	255 x 255 x 30 mm R 393.00	
Silicon mat 	250 x 235 mm R 851.00	Tension clips for towels 	255 x 255 mm R 395.00	

Container complete



590 x 370 x 120 mm
590 x 370 x 180 mm

590 x 370 x 120 mm
590 x 370 x 180 mm

590 x 370 x 120 mm
590 x 370 x 180 mm

590 x 370 x 120 mm
590 x 370 x 180 mm

590 x 370 x 120 mm
590 x 370 x 180 mm

590 x 370 x 120 mm
590 x 370 x 180 mm

Perforated lid
Unperforated bottom

R 110.62
R 110.68
R 120.62
R 120.68
R 130.62
R 130.68
R 140.62
R 140.68
R 150.62
R 150.68
R 160.62
R 160.68

Perforated lid
Perforated bottom

R 111.62
R 111.68
R 121.62
R 121.68
R 131.62
R 131.68
R 141.62
R 141.68
R 151.62
R 151.68
R 161.62
R 161.68

590 x 370 mm

Perforated lid



R 100.51
R 100.52
R 100.53
R 100.54
R 100.55
R 100.56

Unperforated bottom



R 105.62
R 105.68

Perforated bottom



590 x 370 x 120 mm
590 x 370 x 180 mm

R 106.62
R 106.68

Wire baskets



533 x 337 x 50 mm R 190.55


Wire baskets



533 x 337 x 100 mm R 191.60

Container models and dimensions

Perforated lid
Unperforated bottom

300 x 140 x 40 mm  R 410.04 R 420.04 R 430.04 R 440.04 R 450.04 R 460.04		300 x 140 x 70 mm  R 410.07 R 420.07 R 430.07 R 440.07 R 450.07 R 460.07	300 x 140 x 100 mm  R 410.10 R 420.10 R 430.10 R 440.10 R 450.10 R 460.10
300 x 140 x 40 mm Perforated lid Perforated bottom  R 411.04 R 421.04 R 431.04 R 441.04 R 451.04 R 461.04	300 x 140 x 70 mm  R 411.07 R 421.07 R 431.07 R 441.07 R 451.07 R 461.07		
300 x 140 mm Perforated lid  R 400.01 R 400.02 R 400.03 R 400.04 R 400.05 R 400.06	Unperforated bottom  300 x 140 x 40 mm R 405.04 300 x 140 x 70 mm R 405.07	Perforated bottom  300 x 140 x 40 mm R 406.04 300 x 140 x 70 mm R 406.07	
Silicon mat  275 x 125 mm R 850.00	Tension clips for towels  270 x 125 mm R 495.00	Drape  60 x 40 cm R 095.01	

Perforated lid
Unperforated bottom



310 x 190 x 40 mm **R 510.04**

Perforated lid



310 x 190 x 40 mm **R 500.01**

Perforated lid
Unperforated bottom



310 x 190 x 65 mm **R 510.06**

Unperforated bottom



310 x 190 x 40 mm **R 505.04**
310 x 190 x 65 mm **R 505.06**
310 x 190 x 130 mm **R 505.13**

Perforated lid
Unperforated bottom



310 x 190 x 130 mm **R 510.13**

Silicon mat




260 x 180 mm **R 850.10**

Tension clips for towels




277 x 175 mm **R 595.00**

Drape



60 x 40 cm **R 095.01**
65 x 65 cm **R 095.02**
85 x 85 cm **R 095.03**

Optic Container without filter	Perforated lid Unperforated bottom 		Perforated lid Unperforated bottom 
250 x 60 x 50 mm 		450 x 70 x 70 mm 	
Implantat container without filter	Perforated lid Perforated bottom		Perforated lid Perforated bottom
300 x 138 x 65 mm 		500 x 155 x 75 mm 	

Disposable sterile paper filter



R 800.00 Unit 1000 pcs.
for containers 1/1, 3/4, 1/2



R 804.00 Unit 1000 pcs.
for small containers (MINI)



R 802.00 Unit 1000 pcs.
for DENTAL containers

Permanent Textile Filters 50 sterilizations



R 800.01 1 Single
for containers 1/1, 3/4, 1/2




R 804.01 1 Single
for small containers (MINI)




R 802.01 1 Single
for DENTAL containers


Teflon filter 1200 sterilizations



R 800.02 1 Single
for containers 1/1, 3/4, 1/2




R 800.20 1 Single
for small containers (MINI)



R 800.10 1 Single
for DENTAL containers

Paper indicators

69 x 38 mm



R 810.00 100 pcs.
for containers 1/1, 3/4, 1/2

Paper indicators


58 x 18 mm



R 812.00 100 pcs.
for small containers (MINI)

Aluminium indicators without graving

50 x 18 mm



CHIRURGIE	TRANSPLANTATION	GYNÄKOLOGIE
GYNÄKOLOGIE	CHIRURGIE	TRANSPLANTATION
GRUNDSIEB	GRUNDSIEB	GYNÄKOLOGIE
GRUNDSIEB	GRUNDSIEB	GYNÄKOLOGIE

Plastic seals



R 818.00 1000 pcs.




R 830.00
Fixing bolt for divisions

R 830.01
Fixation shaft for instruments



Divisions

R 831.05 50 x 20 mm
R 831.13 130 x 20 mm
R 831.22 225 x 20 mm
R 831.46 460 x 20 mm



R 840.00
Sliders

R 841.00
Rail for sliders

R 842.00
Angle stopper



R 845.00
Spiral spring




R 847.00
Fixation for spiral spring



R 849.00
Fixation clips




Stainless steel perforated sheet baskets



• Stainless steel perforated sheet baskets

89421-50	mm	210	145	50
89424-41	mm	240	240	41
89424-45	mm	240	240	45
89424-50	mm	240	240	50
89424-70	mm	240	240	70
89424-90	mm	240	240	90
89428-40	mm	280	110	40
89428-70	mm	280	110	70
89430-35	mm	300	240	35
89430-50	mm	300	240	50
89430-70	mm	300	240	70
89440-41	mm	380	240	41
89440-45	mm	380	240	45
89440-50	mm	380	240	50
89440-70	mm	380	240	70
89440-90	mm	380	240	90
89448-41	mm	480	240	41
89448-45	mm	480	240	45
89448-50	mm	480	240	50
89448-70	mm	480	240	70
89448-90	mm	480	240	90
89454-41	mm	540	240	41
89454-45	mm	540	240	45
89454-50	mm	540	240	50
89454-70	mm	540	240	70
89454-90	mm	540	240	90
89465-70	mm	650	240	70



• Stainless steel covers for baskets made of perforated stainless steel sheet

89421-14	mm	210	145
89424-24	mm	240	240
89428-11	mm	280	240
89430-24	mm	300	240
89440-24	mm	380	240
89448-24	mm	480	240
89454-24	mm	540	240
89465-24	mm	650	240

Sterilizing cases



- Cover and bottom perforated, with long-life textile filters for steam sterilization

89119-01	220 x 165 x 55 mm
89119-02	300 x 220 x 55 mm
89119-03	400 x 290 x 55 mm




Long-life textile filter



89119-10	178 x 115 mm
89119-20	245 x 165 mm
89119-30	358 x 235 mm

89119-04	220 x 165 x 75 mm
89119-05	300 x 220 x 75 mm
89119-06	400 x 290 x 75 mm

Wire basket for endoscopes



89256-29

290 x 80 x 52 mm
for 1 arthroscope or 1 sinuscope

89257-29

290 x 80 x 52 mm
for 2 arthroscopes or 2 sinuscope

89256-46

460 x 80 x 52 mm
for 1 laparoscope or 1 hysteroscope or 1 cystoscope

89257-46


460 x 80 x 52 mm
for 2 laparoscopes or 2 hysteroscopes or 2 cystoscopes

89256-67

670 x 80 x 52 mm
for 1 bronchoscope

89256-01 / 89256-02 / 89256-03 / 89256-04
Spare parts

Container for Endoscopes



R019-9312

282 x 76 x 33 mm
for 1 arthroscope or 1 sinuscope

R019-9317

432 x 76 x 33 mm
for 1 hysteroscope or 1 cystoscope

R019-9321

511 x 76 x 33 mm
for 1 laparoscope


R019-9326

666 x 114 x 79 mm
for 1 bronchoscope

R019-9331

813 x 132 x 79 mm
for 1 extra long endoscope

Storage- and sterilization rack



89258-47

470 x 250 x 125 mm for 15 MIS instruments
5 / 10 mm Ø

89258-25

232 x 210 x 155 mm for arthroscopy instruments

Silicone bar

Scope Holder



4 cm / 1,5"

R019-5350

15 mm Ø

R019-5353

12 mm Ø



R019-5351

8 mm Ø

R019-5352

12 mm Ø

R019-5354

6 mm Ø

R019-5355

4 mm Ø



7 cm / 2,8"

R019-5340

15 mm/8 mm Ø

R019-5343

15 mm/12 mm Ø



R019-5341

12 mm Ø

Instrument Bar

18 cm / 7.0" + 7 cm / 2,8"



R019-5373 long

R019-5333 short



R019-5374 long

R019-5334 short



R019-5375 long

R019-5335 short

Holder



R019-5210

R019-5200