



HAYDEN MEDICAL



LIPOPLASTY



Company Message



TOM OCASEK
Logistics Manager



HAYDEN CLELLAND
VP Sales & Marketing



SCOTT CLELLAND
Founder, CEO

Hayden Medical is the source for high quality German surgical instruments at reasonable prices. Why buy from the last supplier who charged your hospital, office or clinic a premium, when you can get your best deal right here with Hayden Medical.


Save money today by purchasing from Hayden Medical for all of your specialty plastic and cosmetic surgery instrument needs.

We can personalize your instruments with anything from plain text to the surgery center's logo or the surgeon's signature with our laser marking system. With over 50 years of combined surgical instrument experience, we know which instruments can fit your needs for any cosmetic or plastic surgery procedure.

We will chase down your hard-to-find items and take the guesswork out of your day. High quality German instrumentation, great service and affordable pricing is our specialty.

If you've used Hayden instruments for years or are a first time customer, you will see and feel the difference every time one of our instruments are in your hands.

RETURNED GOODS POLICY



All returns must be accompanied by a Returned Goods Authorization number, available by calling our Customer Service Department.

THE FOLLOWING INFORMATION IS NEEDED TO REQUEST AUTHORIZATION FOR A RETURN.

- REASON FOR RETURN
- DATE OF PURCHASE
- INVOICE OR PO NUMBER
- PRODUCT NUMBER AND SIZE
- QUANTITY TO RETURN

All authorized returns must be made within thirty (30) days of the invoice date and accompanied by the appropriate purchase order number and Returned Goods Authorization number.

ALLOWABLE SCHEDULE FOR AUTHORIZED RETURN.

- Products processed in error by Hayden Medical, Inc. will be credited in full
- Defective product will be credited in full
- Unless defective or shipped in error, any return accepted is subject to up to a 25% restocking fee (minimum of \$10.00)
- All returned items are subject to Hayden Medical, Inc. inspection and acceptance
- Product must be shipped prepaid

NON-RETURNABLE, OTHER THAN FOR MANUFACTURING DEFECT.

- Engraved instruments or instruments etched or labeled by the hospital, surgery center or surgeon's office
- Products damaged beyond repair or refurbishment or not in saleable condition
- Instruments returned more than thirty (30) days after invoice date
- Custom, laser marked, special or modified items
- Disposable items and Cannulas

Used products must be returned sterile wrapped or they will not be accepted.



SEND ALL RETURNS TO:
19425 Soledad Cyn. Road, #411
Santa Clarita, CA 91351



CANNULA SETS

LIPOPLASTY



Micro Aspiration Harvesting
Style 1 Micro Injector
Style 2 Micro Injector
Style 3 Micro Injector
Facial Anesthesia Needle
Snapper Lock for 20cc Syringe
Micro Transfer Hub Set

AUTOLOGOUS FAT TRANSFER SET (FACE)

ORDER ENTIRE SET
HC-FACE

| PRODUCT CODE | DESCRIPTION |
|--------------|---|
| HC-MAN1215 | Micro Aspiration Harvesting Cannula - 12ga x 15cm |
| HC-MIN11607 | Style 1 Micro-Injector - 16ga x 7cm |
| HC-MIN21609 | Style 2 Micro-Injector - 16ga x 9cm |
| HC-MIN31609 | Style 3 Micro-Injector - 16ga x 9cm |
| HC-FAN1615 | Facial Anesthesia Needle - 16ga x 15cm |
| HC-SL20 | Snapper Lock for 20cc Syringe |
| NFT-SET | Micro Transfer Hub Set - ID 1.2mm, 1.4mm & 2.4mm |



Tumescence Infiltration Cannula
Tri-Port Harvester
Style 1 Micro Injector
Snapper Lock for 60cc Syringe
Transfer Hub Female Luer Ends

AUTOLOGOUS FAT TRANSFER SET (BREAST)

ORDER ENTIRE SET
HC-BRST

| PRODUCT CODE | DESCRIPTION |
|--------------|---|
| HC-BAN1215 | Tumescence Infiltration Cannula - 12ga x 15cm |
| HC-III0315 | Tri-Port Harvester - 3mm x 15cm |
| HC-MIN11215 | Style 1 Micro-Injector - 12ga x 15cm |
| HC-SL60 | Snapper Lock for 60cc Syringe |
| NFT-2.4 | Transfer Hub, Female Luer Ends - ID 2.4mm |



Mercedes Cannula (2)
Style 1 Injector (2)
Tumescence Infiltration Cannula

AUTOLOGOUS FAT TRANSFER SET (GLUTEAL)

ORDER ENTIRE SET
HC-GLUT

| PRODUCT CODE | DESCRIPTION |
|--------------|---|
| HC-MER0420 | Mercedes Cannula - 4mm x 20cm |
| HC-MER0425 | Mercedes Cannula - 4mm x 25cm |
| HC-MIN10315 | Style 1 Injector - 3mm x 15cm |
| HC-MIN10320 | Style 1 Injector - 3mm x 20cm |
| HC-BAN1225 | Tumescence Infiltration Cannula - 12ga x 25cm |



Tumescence Infiltration Cannula (2)

Black Infusion Handle

TUMESCENT INFILTRATION SET

ORDER ENTIRE SET
HC-TINF

| PRODUCT CODE | DESCRIPTION |
|---------------|---|
| HC-BAN1232 | Tumescence Infiltration Cannula - 12ga x 32cm |
| HC-BAN1415 | Tumescence Infiltration Cannula - 14ga x 15cm |
| HC-INFU-C-BLK | Infusion Handle with On/Off Control - Black |

LIPOPLASTY

Hayden Medical's cannulas and cannula accessories are made of high-quality, medical-grade stainless steel and made to order in the USA. If there is a size or pattern that you don't see, we can make it. Please contact us directly with custom cannula requests.



HAYDEN MULTI-HOLE HARVESTER

270° of extraction and 18 holes for low-volume liposuction procedures (available with angled shaft). Typically used to harvest for small-volume fat grafting to areas like the face and hands.

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|--------------|
| HC-SAT0315 | 3mm | 15cm | LL / LH / TM |
| HC-SAT0326 | 3mm | 26cm | LL / LH / TM |
| HC-SAT0332 | 3mm | 32cm | LL / LH / TM |



MERCEDES

360° of extraction and 3 holes for high-volume liposuction procedures. Typically the most popular and commonly used liposuction cannula.

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|--------------|
| HC-MER0326 | 3mm | 26cm | LL / LH / TM |
| HC-MER0332 | 3mm | 32cm | LL / LH / TM |
| HC-MER0426 | 4mm | 26cm | LL / LH / TM |
| HC-MER0432 | 4mm | 32cm | LL / LH / TM |
| HC-MER0526 | 5mm | 26cm | LL / LH / TM |
| HC-MER0532 | 5mm | 32cm | LL / LH / TM |



BASKET

360° of extraction and 3 flared openings for large-volume liposuction procedures. These can be referred to as a "Flared Mercedes" cannula as well.

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|--------------|
| HC-BAS0326 | 3mm | 26cm | LL / LH / TM |
| HC-BAS0332 | 3mm | 32cm | LL / LH / TM |
| HC-BAS0426 | 4mm | 26cm | LL / LH / TM |
| HC-BAS0432 | 4mm | 32cm | LL / LH / TM |
| HC-BAS0526 | 5mm | 26cm | LL / LH / TM |
| HC-BAS0532 | 5mm | 32cm | LL / LH / TM |



MERCEDES PLUS (DOUBLE MERCEDES)

360° of extraction and 6 holes for high volume liposuction procedures.

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|--------------|
| HC-MER+0326 | 3mm | 26cm | LL / LH / TM |
| HC-MER+0332 | 3mm | 32cm | LL / LH / TM |
| HC-MER+0426 | 4mm | 26cm | LL / LH / TM |
| HC-MER+0432 | 4mm | 32cm | LL / LH / TM |
| HC-MER+0526 | 5mm | 26cm | LL / LH / TM |
| HC-MER+0532 | 5mm | 32cm | LL / LH / TM |



CANDY CANE

360° of extraction for medium to high-volume liposuction procedures.

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|--------------|
| HC-CC0326 | 3mm | 26cm | LL / LH / TM |
| HC-CC0332 | 3mm | 32cm | LL / LH / TM |
| HC-CC0426 | 4mm | 26cm | LL / LH / TM |
| HC-CC0432 | 4mm | 32cm | LL / LH / TM |
| HC-CC0526 | 5mm | 26cm | LL / LH / TM |
| HC-CC0532 | 5mm | 32cm | LL / LH / TM |

LE PATTERNS

LIPOPLASTY




SPATULA

90° of extraction and a single hole for low to medium-volume liposuction procedures. Non-traumatic spatula tips are ideal for delicate areas like the face and neck.

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|--------------|
| HC-SPA0215 | 2mm | 15cm | LL / LH / TM |
| HC-SPA0315 | 3mm | 15cm | LL / LH / TM |
| HC-SPA0415 | 4mm | 15cm | LL / LH / TM |
| HC-SPA0515 | 5mm | 15cm | LL / LH / TM |

HUB OPTIONS



SPATULA 2-PORT

90° of extraction and 2 holes for low to medium-volume liposuction procedures. Non-traumatic spatula tips are ideal for delicate areas like the face and neck.

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|--------------|
| HC-SP20215 | 2mm | 15cm | LL / LH / TM |
| HC-SP20315 | 3mm | 15cm | LL / LH / TM |
| HC-SP20415 | 4mm | 15cm | LL / LH / TM |
| HC-SP20515 | 5mm | 15cm | LL / LH / TM |



MODERN SPIRAL

270° of extraction with multiple holes for high-volume aspiration. The Modern Spiral cannula is becoming the new favorite of surgeons who traditionally use a Tri-Port or Mercedes cannula.

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|--------------|
| HC-MS0326 | 3mm | 26cm | LL / LH / TM |
| HC-MS0332 | 3mm | 32cm | LL / LH / TM |
| HC-MS0426 | 4mm | 26cm | LL / LH / TM |
| HC-MS0432 | 4mm | 32cm | LL / LH / TM |
| HC-MS0526 | 5mm | 26cm | LL / LH / TM |
| HC-MS0532 | 5mm | 32cm | LL / LH / TM |



TRI-PORT

270° of extraction and 3 holes for high-volume liposuction procedures

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|--------------|
| HC-III0215 | 2mm | 15cm | LL / LH / TM |
| HC-III0226 | 2mm | 26cm | LL / LH / TM |
| HC-III0315 | 3mm | 15cm | LL / LH / TM |
| HC-III0326 | 3mm | 26cm | LL / LH / TM |
| HC-III0332 | 3mm | 32cm | LL / LH / TM |
| HC-III0415 | 4mm | 15cm | LL / LH / TM |
| HC-III0426 | 4mm | 26cm | LL / LH / TM |
| HC-III0432 | 4mm | 32cm | LL / LH / TM |
| HC-III0526 | 5mm | 26cm | LL / LH / TM |
| HC-III0532 | 5mm | 32cm | LL / LH / TM |



TRI-PORT PLUS FOR LARGE-VOLUME HARVESTING

270° of extraction and 6 holes for high-volume liposuction procedures.

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|--------------|
| HC-III+0326 | 3mm | 26cm | LL / LH / TM |
| HC-III+0332 | 3mm | 32cm | LL / LH / TM |
| HC-III+0426 | 4mm | 26cm | LL / LH / TM |
| HC-III+0432 | 4mm | 32cm | LL / LH / TM |
| HC-III+0526 | 5mm | 26cm | LL / LH / TM |
| HC-III+0532 | 5mm | 32cm | LL / LH / TM |



INJECTION & INFILT

MICRO INJECTION



STYLE 1

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|-------------|
| HC-MIN11607 | 16ga | 7cm | LUER LOCK |
| HC-MIN11609 | 16ga | 9cm | LUER LOCK |
| HC-MIN11807 | 18ga | 7cm | LUER LOCK |
| HC-MIN11809 | 18ga | 9cm | LUER LOCK |
| HC-MIN11905 | 19ga | 5cm | LUER LOCK |
| HC-MIN12003 | 20ga | 3cm | LUER LOCK |
| HC-MIN12005 | 20ga | 5cm | LUER LOCK |



STYLE 2

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|-------------|
| HC-MIN21607 | 16ga | 7cm | LUER LOCK |
| HC-MIN21609 | 16ga | 9cm | LUER LOCK |
| HC-MIN21807 | 18ga | 7cm | LUER LOCK |
| HC-MIN21809 | 18ga | 9cm | LUER LOCK |



STYLE 3

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|-------------|
| HC-MIN31607 | 16ga | 7cm | LUER LOCK |
| HC-MIN31609 | 16ga | 9cm | LUER LOCK |
| HC-MIN31807 | 18ga | 7cm | LUER LOCK |
| HC-MIN31809 | 18ga | 9cm | LUER LOCK |



"V" DISSECTOR

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|-------------|
| HC-MINV1215 | 12ga | 15cm | LUER LOCK |
| HC-MINV1225 | 12ga | 25cm | LUER LOCK |
| HC-MINV1415 | 14ga | 15cm | LUER LOCK |
| HC-MINV1409 | 14ga | 9cm | LUER LOCK |
| HC-MINV1607 | 16ga | 7cm | LUER LOCK |
| HC-MINV1609 | 16ga | 9cm | LUER LOCK |
| HC-MINV1807 | 18ga | 7cm | LUER LOCK |

MICRO ASPIRATION



MERCEDES

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|-------------|
| HC-MER1215 | 12ga | 15cm | LUER LOCK |
| HC-MER1225 | 12ga | 25cm | LUER LOCK |
| HC-MER1415 | 14ga | 15cm | LUER LOCK |
| HC-MER1425 | 14ga | 25cm | LUER LOCK |
| HC-MER0315 | 3mm | 15cm | LUER LOCK |
| HC-MER0325 | 3mm | 25cm | LUER LOCK |



FAT ASPIRATION CANNULA (FAC)

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|-------------|
| HC-MAN1215 | 12ga | 15cm | LUER LOCK |
| HC-MAN1220 | 12ga | 20cm | LUER LOCK |
| HC-MAN1225 | 12ga | 25cm | LUER LOCK |
| HC-MAN1415 | 14ga | 15cm | LUER LOCK |
| HC-MAN1425 | 14ga | 25cm | LUER LOCK |

RATION CANNULAS

BODY TISSUE INJECTION

TUMESCENT INFILTRATION



STYLE 1

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|-------------|
| HC-MIN11215 | 12ga | 15cm | LUER LOCK |
| HC-MIN11220 | 12ga | 20cm | LUER LOCK |
| HC-MIN11225 | 12ga | 25cm | LUER LOCK |
| HC-MIN11415 | 14ga | 15cm | LUER LOCK |
| HC-MIN11420 | 14ga | 20cm | LUER LOCK |
| HC-MIN11425 | 14ga | 25cm | LUER LOCK |




BODY ANESTHESIA NEEDLE

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|-------------|
| HC-BAN1015 | 10ga | 15cm | LUER LOCK |
| HC-BAN1032 | 10ga | 32cm | LUER LOCK |
| HC-BAN1215 | 12ga | 15cm | LUER LOCK |
| HC-BAN1225 | 12ga | 25cm | LUER LOCK |
| HC-BAN1232 | 12ga | 32cm | LUER LOCK |
| HC-BAN1415 | 14ga | 15cm | LUER LOCK |
| HC-BAN1425 | 14ga | 25cm | LUER LOCK |
| HC-BAN1432 | 14ga | 32cm | LUER LOCK |
| HC-BAN1615 | 16ga | 15cm | LUER LOCK |

TISSUE TRANSFER

SOFT TISSUE TRANSFER SYSTEM

SYRINGE, LUER LOCK TO LUER LOCK



| PRODUCT CODE | DIAMETER |
|--------------|----------|
| HC-LL3 | 3mm |
| HC-LL4 | 4mm |



FACIAL ANESTHESIA NEEDLE

| PRODUCT CODE | DIAMETER | LENGTH | HUB OPTIONS |
|--------------|----------|--------|-------------|
| HC-FAN1415 | 14ga | 15cm | LUER LOCK |
| HC-FAN1615 | 16ga | 15cm | LUER LOCK |

MICRO TRANSFER HUB SET



1.2mm



1.4mm



2.4mm

| PRODUCT CODE | DESCRIPTION |
|--------------|-------------|
| NFT-SET | Set of 3 |




TWIST LOCK

| PRODUCT CODE | DESCRIPTION |
|--------------|---|
| HC-MJ-TL60 | for 60cc Monoject Syringe |
| HC-BD-TL60 | for 60cc BD Syringe |
| HC-MJ-TL35 | for 35cc Monoject Syringe |
| HC-BT-TL20 | for 20cc BD and 20cc Terumo Syringe |
| HC-BD-TL30 | for 30cc BD Syringe |
| HC-BD-TL10 | for 10cc BD Syringe |
| HC-MJ-TL10 | for 10cc Monoject and 12cc Monoject Syringe |



LIPOSUCTION HANDLES & ACCESSORIES

INFUSION HANDLE 3/4" with On/Off Control



| PRODUCT CODE | DESCRIPTION |
|---------------|---|
| HC-INFU-C-BLK | A flow on/off control infusion handle for use with roller pump devices. |


INFUSION HANDLE



| PRODUCT CODE | DESCRIPTION |
|--------------|--|
| HC-INFU-NC-L | Standard Infusion handle, No Control, Ideal for infusion needle 1" diameter. |

PLASTIC STERILE TRAYS

- Autoclavable
- Silicone Finger-Tip Mats



| PRODUCT CODE | DESCRIPTION |
|--------------|---------------------------------|
| 030-105 | 4" x 7-1/2" x 3/4" |
| 030-105.5 | 4" x 7-1/2" x 1-1/2" |
| 030-107 | 6" x 10" x 3/4" |
| 030-108 | 6" x 10" x 1-1/2" |
| 030-108D | 6" x 10" x 1-1/2" Double-Decker |
| 030-111 | 10" x 15" x 3/4" |
| 030-112 | 10" x 15" x 1-1/2" |

NEEDLE CONNECTOR



| PRODUCT CODE | DESCRIPTION |
|--------------|--------------------------------------|
| HC-NC | Needle connector for use with tubing |

ASPIRATION HANDLES



| PRODUCT CODE | DESCRIPTION |
|--------------|---------------------------------|
| HC-AH-S | Standard Aspiration Handle 3/4" |
| HC-AH-L | Standard Aspiration Handle 1" |

AWLS



| PRODUCT CODE | DESCRIPTION |
|--------------|-------------|
| HC-EP3 | 3mm |

SNAPPER LOCKS



| PRODUCT CODE | DESCRIPTION |
|--------------|---|
| HC-BT-SL10 | for 10cc BD and 10cc Terumo Syringe |
| HC-MJ-SL10 | for 10cc Monoject and 12cc Monoject Syringe |
| HC-BD-SL20 | for 20cc BD Syringe |
| HC-MJ-SL20 | for 20cc Monoject Syringe |
| HC-BD-SL30 | for 30cc BD Syringe |
| HC-MJ-SL35 | for 35cc Monoject Syringe |
| HC-TR-SL60 | for 60cc Terumo Syringe |
| HC-MB-SL60 | for 60cc Monoject and B-D Syringe |

ADIPOSE TISSUE PROCESSING

PRISM PARTICLE SIZER™

The **PRISM PARTICLE SIZER™** is a particle-sizing and filtering device for human tissue and/or fluid. When used to process fat aspirate, the device creates precision-sized micro-fragmented adipose tissue (MFAT) across a selectable dynamic sizing range of 3000-µm to 200-µm.

BENEFITS OF PPS DEVICE INCLUDES:

- More ergonomic single pass sizing & filtering
- Creates the most purified processed fat aspirate
- Cost-effective, reusable, Class I medical device
- Patent pending and manufactured in the USA
- Recommended PPS device accessories include:

PRISM CAROUSEL™

PRISM TRANSFER CANNULA™

PRISM CANNULA CLEANER™



| PRODUCT CODE | DESCRIPTION | | | | | |
|------------------|-------------|---------|---------|---------|---------|--|
| PPS-Face | 2000-µm | 1500-µm | 900-µm | 500-µm | 300-µm | |
| PPS-Breast | 2500-µm | 1500-µm | 1200-µm | | | |
| PPS-Hands | 2000-µm | 1200-µm | 1000-µm | | | |
| PPS-Buttocks | 3000-µm | 2500-µm | 2000-µm | 1500-µm | | |
| PPS-Regenerative | 3000-µm | 2500-µm | 2000-µm | 1500-µm | | |
| PPS-Set | 3000-µm | 2500-µm | 2000-µm | 1500-µm | 1200-µm | |
| | 1000-µm | 900-µm | 800-µm | 700-µm | 600-µm | |
| | 500-µm | 400-µm | 300-µm | 200-µm | | |

PRISM PARTICLE SIZER ACCESSORIES

PRISM CAROUSEL™




| PRODUCT CODE | DESCRIPTION |
|--------------|--|
| P-CAR-6 | One Decanting/Storage Rack for 6 Luer-Lock Syringes |
| P-CAR-12 | One Decanting/Storage Rack for 12 Luer-Lock Syringes |

PRISM CANNULA CLEANER™



| PRODUCT CODE | DESCRIPTION |
|--------------|---|
| PCC | One reusable nylon tool for cleaning and removing debris when reprocessing PPS device |

PRISM TRANSFER CANNULA™



| PRODUCT CODE | DESCRIPTION |
|--------------|---|
| PTC-SET | Set of one shorter and one longer Cannula for transferring fat aspirate while using the PPS device system |

The **PRISM PARTICLE SIZER™** medical device and its recommended accessories are FDA-registered and manufactured by HOGUE SURGICAL (Maple Grove, MN, USA)

PUMPS & CANISTER

INFILTRATION PUMP

PRODUCT CODE

HC-IP

- Delivers High-Volume of Tumescent Fluid
- Can Sit on Aspirator
- Compact and Minimal Noise



ASPIRATOR PUMP

PRODUCT CODE

HC-AP4

- High-Torque Suction
- Portable and Compact
- Integrated Pressure Gauge
- Minimal Device Noise



FAT COLLECTION CANISTER KIT

- Autoclavable canisters for use during liposuction procedures where fluid/tissue needs to be extracted from the bottom of the canister. Extraction is performed with a luer lock syringe by connecting it to the reusable luer lock extraction tube provided in the complete kit.

Canister Kits are Completely autoclavable and include:

- Fat Transfer Canister
- Canister Lid
- Extraction Tube
- Canister Stand
- Tube Adapter

| PRODUCT CODE | DESCRIPTION |
|--------------|-------------|
| HFAT-1000CC | 1000cc Kit |
| HFAT-3000CC | 3000cc Kit |



LASER & FACIAL PLASTIC SURGERY

OCULAR LASER SHIELD WITH SUCTION CUP



| PRODUCT CODE | DESCRIPTION | | | |
|--------------|-------------|---------------|------|--|
| OP21-261P | X-Small | 22.0 x 19.5mm | Pair | |
| OP21-262P | Small | 24.0 x 21.0mm | Pair | |
| OP21-264P | Medium | 25.5 x 22.0mm | Pair | |
| OP21-266P | Large | 26.5 x 23.0mm | Pair | |

OCULAR LASER SHIELD WITH OFFSET HANDLE



| PRODUCT CODE | DESCRIPTION | | | |
|--------------|-------------|---------------|------|--|
| OP21-261HP | X-Small | 22.0 x 19.5mm | Pair | |
| OP21-262HP | Small | 24.0 x 21.0mm | Pair | |
| OP21-264HP | Medium | 25.5 x 22.0mm | Pair | |
| OP21-266HP | Large | 26.5 x 23.0mm | Pair | |

PLASTIC OCULAR SHIELD WITH SUCTION CUP



| PRODUCT CODE | DESCRIPTION | | | |
|--------------|-------------|---------------|------|--|
| OP23-201BP | X-Small | 23.5 x 21.5mm | Pair | |
| OP23-202BP | Small | 26.0 x 23.5mm | Pair | |
| OP23-204BP | Medium | 27.5 x 24.5mm | Pair | |
| OP23-206BP | Large | 28.5 x 25.5mm | Pair | |

PLASTIC OCULAR SHIELD WITH HANDLE



| PRODUCT CODE | DESCRIPTION | | | |
|--------------|-------------|---------------|------|--|
| OP23-201BHP | X-Small | 23.5 x 21.5mm | Pair | |
| OP23-202BHP | Small | 26.0 x 23.5mm | Pair | |
| OP23-204BHP | Medium | 27.5 x 24.5mm | Pair | |
| OP23-206BHP | Large | 28.5 x 25.5mm | Pair | |

SUBMENTAL FIBEROPTIC RETRACTOR



| PRODUCT CODE | DESCRIPTION |
|--------------|--|
| 210-180L | 4-1/4" , 49mm Tapered Blade, 60mm x 18mm |

DESMARRES RETRACTOR NON-CONDUCTIVE




| PRODUCT CODE | DESCRIPTION |
|--------------|------------------------------|
| 210-194NC | 11mm wide |
| 210-195NC | 13mm wide |
| 210-196NC | 15mm wide |
| 210-197NC | 17mm wide |
| 210-198NC | 8 / 10mm wide - double ended |

PAR SCISSORS SUPER-CUT



| PRODUCT CODE | DESCRIPTION |
|--------------|--------------------------------|
| 100-321SC | Curved, serrated blade, 4-1/2" |

ADSON FORCEPS, WITH TUNGSTEN INSERT



| PRODUCT CODE | DESCRIPTION |
|--------------|---|
| 540-136 | Tissue, 1x2 teeth, serrated, 1.4mm standard |
| 540-139F | Tissue, 1x2 teeth, serrated, 1mm delicate |
| 540-117 | Dressing, serrated, 1.4mm standard |
| 540-124 | Brown tip, serrated, 3x3 teeth |



SURGICAL SUITE EQUIPMENT

Outfit your surgical suite and procedure rooms with our wide range of new and refurbished equipment. We offer major brands and can match your needs with multiple choices to suit your specific budget. We pride ourselves on the quality of our products, excellent customer service and service after the sale.

• SURGICAL LIGHTING

0ALED070SC



0ALED0100DC



• STAINLESS STEEL TABLES AND STANDS

HCM731



HCM513



HCM520



• PROCEDURE CHAIRS

ESC2



• PATIENT MONITORS



• TABLE TOP AUTOCLAVES



• DOCTOR'S STOOL

EDS-1060G



• SURGICAL STOOL

S900-W



ANESTHESIA MACHINES, AUTOCLAVES, SURGICAL C-ARMS, C-ARM TABLES, DEFIBRILLATORS, ELECTROSURGICAL UNITS, ENDOSCOPY EQUIPMENT, EXAM TABLES, MEDICAL STAINLESS STEEL, MEDICAL STRETCHERS, PATIENT MONITORS, PROCEDURE TABLES, SURGICAL MICROSCOPES, SURGICAL LIGHTING, SURGICAL TABLES, WARMING CABINETS



SmoothWave i200

4.0MHz RF technology with Bipolar Turbo and Tissue Matching Intelligence provides greater surgical precision and soft tissue preservation.



3-Button Fingerswitch Wand Disposable (10pcs)
SIEC-3FHPB/D



Microincision Tip, 5mm Needle Tip Length,
3cm Shaft Length (2pcs)
STEE301-2



Fine Wire, 3cm Shaft Length (25pcs) Hex Shaft
SST-W-25



Broad Needle (25pcs) Hex Shaft
SBN-25

